# BEAUTIFLY

Meso Intense

# Golden

series

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Power supply:	AC/DC
Input:	AC 100-240V 50/60Hz
Output:	DC-5V, 0.5A
Rotate speed:	10000-15000 RPM
Needle length:	0-2.5MM
Speed adjustable:	6 levels
Included:	Device, charger, charging cable, user manual, 3 pcs of 12-needle cartridges

### Security

Before using the product, read this manual and follow the instructions contained therein during use of devices.

Device description [ see page 2, f1 ]

1. Power button
2. Photon light switch
3. Mode display
4. Needle length adjustment ring
5. Charging input

### Needle cartridge installation

1. Remove the needle cartridge from the sterilized package.
2. Match round slot 1 to 2 of the machine and push the needle cartridge in horizontally as far as it will go.
3. Turn the needle cartridge clockwise to lock it.
4. To remove the needle cartridge, turn it clockwise and push it out horizontally.

### Needle depth adjustment

Adjust the depth of the needles by turning the scale adjustment ring. Turn it clockwise to shorten the needles and counterclockwise to lengthen them.

### Instructions for use

1. Wash and dry your hands thoroughly before performing the procedure.
2. Use a gentle cleanser to thoroughly cleanse your face.
3. Apply the appropriate ampoule or preparation for microneedle mesotherapy.
4. Remove the needle cartridge from the sterile package, attach it to the device and adjust the needle depth (recommended needle length for first use is 0.25mm).
5. Remove the protective cap from the needle cartridge and dispose of it in general household waste.
6. Press and hold the power button for 1.5 seconds to start up the device.
7. Short press the button to switch gears. Your run information will be displayed on the screen.
8. Short press the phototherapy switch to select one of the 4 lights (green, red, blue, yellow).
9. Perform the treatment for 10-15 minutes, holding the device in a vertical position.
10. Gently massage the excess cosmetics left on the skin after the treatment.

## Charging

Before the first treatment, make sure that the device is fully charged.

Use **only** the USB charger provided with the product.

When power off and charging, the display **will flash** and stop flashing when fully charged.

Do not use the charger near a TV or radio.

## Cleaning and maintenance

Store the device in an upright position to prevent cosmetics from entering it.

Wipe with a damp cloth after use.

Do not immerse the device in water, do not hold it under a stream of water - the device is not waterproof.

## Warnings

- The needle cartridge must be moved continuously and may be stopped on the same area of skin for 3-5 seconds, no more than twice.
- Temporary redness may appear on the skin after using the device. If these symptoms do not improve over time, consult your doctor.
- Avoid skin exposure to the sun, eating spicy foods, seafood, alcoholic beverages and other stimulating foods up to three days after surgery.
- Needle cartridges can **only** be used once and should be disposed of safely after use.
- Needle cartridges are for personal use only. They should not be shared with others to avoid infection.
- Do not use the needle cartridge if the disposable package is damaged or cracked.
- Perform a skin test on your hand before first use. Do not use the device if your skin itches or has an allergic reaction.
- Parts of the device other than the needle tips must not come into contact with the liquid to avoid damage to the device.
- The needle cartridge can **only** be used by one person.
- Do not share the same needle cartridge with other people or reuse it. Otherwise, it may cause skin hypersensitivity, irritation, inflammation and other problems. In the event of a malfunction, do not attempt to repair or disassemble the device yourself.

## Storage

- Store the device without the needle cartridge installed.
- Store in a dry and clean place.
- Store away from high temperatures.
- Keep out of reach of children.
- Place the device in a dry, ventilated place.
- Do not expose the device to direct sunlight.

More information about the product is available on the manufacturer's website at:  
[ see page 2, f2 ]

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EN: Electrical and electronic devices, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out trash can symbol on the packaging. Used electrical and electronic devices, batteries, and accumulators must be taken to local selective waste collection points or to the retailer, and more information can be found in your municipality.

Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components harmful to the environment and human health. Improper handling or damage to them can lead to health hazards or environmental pollution during further disposal or recycling.

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EN: The manufacturer hereby declares, on their sole responsibility, that the product "Meso Intense Golden" complies with the Directive 2014/30/EU of the European Parliament and of the Council of 26th February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of 8th June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, along with Commission Delegated Directive (EU) 2015/863 of 31st March 2015 amending Annex I to Directive 2011/65/EU with regard to the list of restricted substances. The full text of the Declaration of Conformity is available at the website:

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