BEAUTIFLY

Scanner Smart App

Body
series

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- www.beautifyfy.eu
- www.youtube.com/beautifyfyeu

Type of scale:	Electronic
Surface type:	Glass
Maximum weight:	180kg
Power supply:	4xAAA batteries
Functions:	internal memory, 29 measurements, including: body analysis, weight, fat tissue content, body water, muscle measurement, bone mass, BMI, Kcal.

Safety

Before using the product, read the instructions and follow the recommendations contained in it during the use of the device.

Device description [page 2 - f1]

1. LED display
2. USB cable
3. Handheld grip
4. Electrodes

Display description [page 2 - f2]

1. Display data / body mass
2. Bluetooth
3. Units
4. Male/Female/Athlete
5. User profile
6. Fat tissue level (%)
7. Muscles (%)
8. BMI
9. Fat tissue level analysis: low/standard/moderate/overweight/obesity
10. Body water (%)
11. Bone mass
12. Kcal

How to use

1. Before first use, install batteries. Open the cover on the back of the scale and insert the batteries according to the polarity indications (+/-) in the battery compartment.
2. Place the scale on a hard and flat surface. Do NOT use it on a carpet or other soft surface. Make sure that your feet and the surface of the scale are dry to avoid slipping.
3. Stand straight with bare feet, evenly covering the surface of the scale [page 2 - f3]
4. The scale offers two measurement options: using 4 electrodes - position 1 [page 2 - f4]. For a fuller data measurement, hold the grip in your hands (in two places) - position 2 [page 2 - f4]
5. Press the "SET", "▲", or "▼" button on the scale to turn it on, and the scale will move to the user parameter settings. The device stores personal parameters for 9 people (P1-P9).
6. Select your profile.

Note: the scale has an identification function, there is no need to set user data every time. If the weight difference between you and other users is within +/-3 kg during one weighing, the scale will identify all users and display the corresponding user number.
7. Stand still for a few seconds and do not step off the scale until the data is displayed on the LCD. The device displays indicators: weight, fat content, body water, muscles, bone mass, BMI, Kcal.

Note

- The scale will automatically turn off after 20 seconds if no other operation is performed.
- When the scale is not connected to Bluetooth - the LCD will display the user number from the range "P1" - "P9"; if Bluetooth is connected, the LCD will display the user number "P0".

Setting personal parameters

1. press <SET> to turn on the scale.
2. press ▲ or ▼ to switch within P1-P9 and confirm by pressing <SET>
3. to set personal parameters press ▲ or ▼ and select one of the following options, then confirm by pressing <SET>:
 - male - female - athlete,
 - height,
 - age
4. press ▲ or ▼ to select the mass unit (kg/lb/st.lb). Press <SET> to confirm.

Using the application

1. Download the application and register
The scale uses the Ailink application, which can be downloaded on iOS <App Store> and Google Play platform or you can also scan the QR code [page 2 - f5] to download the application. After downloading, register an account.
2. Pair the device

Before pairing, make sure Bluetooth is turned on
Then go to the <me> option to complete your personal data.
3. Turn on the scale and open the Ailink application. Choose the option <Add manually> or <Nearby Devices> to pair devices.
- choose <Nearby Devices>, then touch <Fat Tissue Scale>(Bluetooth)
- choose the <Add manually> option, and then touch <Fat Tissue Scale> (Bluetooth)
To finish the setup, follow the instructions in the app until the message <measurement> appears

Note: during the pairing process, keep the scale turned on.

4. Body composition measurement
Step on the scale, which will automatically synchronize with your device. The scale will perform a body composition analysis and display measurements in the app.

- You can check your measurements in the <History> tab
- In the <Buddy> tab, you can add your dietitian or fitness trainer - with each measurement, they will receive data about your body.
- Depending on the mobile phone, you can synchronize your data with Google Fit, Fitbit, or Apple Health, in the <settings> tab in your profile.
- If you use the scale without turning on the application, the data will be saved and will wait for synchronization during the next application launch.

Troubleshooting

1. When the LED displays the message "LO", install a new battery.
2. When the LED shows "Err", it means that the maximum weight has been exceeded.
3. When the LED shows "Err1", it means a measurement error. To eliminate it, perform one of the following actions:
 - remove socks
 - stand in the correct position
 - check the information about the settings, provide the correct personal data.

Dangers and warnings

1. Do not use replacement parts not approved by the manufacturer.
2. Non-professional service should not disassemble and replace parts of the device without permission.
3. If the scale is malfunctioning, check the batteries. If necessary, replace the batteries.
4. If the percentage of body fat is less than 5% or more than 50% the scale shows an error - check if the user information is correct.
5. The scale is a precise measuring device. Never jump or hit the scale.
6. Store away from water, heat, and extreme cold.
7. Make sure the device is not exposed to direct sunlight or placed near radiators and other heat sources.
8. The data provided on the scale is for reference only. Before starting any diet or exercise program, consult a doctor.
9. The scale is not a toy
10. The device is not suitable for pregnant women and people with a pacemaker

Cleaning and maintenance

Before cleaning, turn off the product. To clean a dirty surface, use a soft cloth moistened with alcohol or glass cleaner. Do not clean the product and accessories with corrosive detergents, such as benzene, thinners. Store in a cool, dry place when not in use. Also, avoid exposure to direct sunlight or a humid environment, otherwise, it may cause fading of colors or damage to the product.

More information about the company and products [page 2 - f6]

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after 13.08.2005.

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[EN] The manufacturer declares under its sole responsibility that the Smart Scanner Smart App product is in compliance with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, Directive of the European Parliament and of the Council 2014/53/EU of 16 April 2014 on the harmonization of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list restricted substances. The full content of the declaration of conformity is available at: www.beautify.eu/deklaracja

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