# BEAUTIFLY

Derma Ice

PRO  
series

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Product Size:	190*45*40mm
Product weight:	155g
Working Voltage:	3.7V
Charging voltage:	5V
Material:	PC+ABS
Suction power:	61Kpa±5Kpa
Battery Capacity:	1000mAh
Charging time:	150min

### Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

### Device description [page 2-f1]:

1. Replaceable head.
2. On/off button intensity adjustment.
3. Hot/ice mode switch button, hot and cold compress.
4. Charging port.
5. Air intake.

### Part description:

1. Microcrystalline head: Microcrystalline tip with diamond particles, removes dead and calloused epidermis, acne scars.
2. Large round head: Vacuum lifting produces a gentle, stimulating massage effect, improves microcirculation, skin becomes tighter.
3. Small round head: Makes it easier to remove blackheads and blackheads from hard-to-reach places.
4. Oval head: Used for fine facial lines, lifts and strengthens the structure and elasticity of the skin.

### Principles of operation:

1. Lifting massage: The vacuum massage system improves the firmness and oval of the face. Vacuum massage supports better blood and lymph circulation, making the skin tighter, better supplied with blood, and reducing swelling and sagging skin. Vacuum absorption technology, used for skin lifting and massage, increases the absorption of cosmetics. This makes the complexion brighter and radiant. At the same time, it strengthens the skin's defense capabilities and keeps it firm and elastic.
2. Microcrystalline head: Produces over 100,000 natural abrasive particles that thoroughly remove dead skin from the facial surface. The skin remains smooth and delicate.
3. Blue photon light: Designed for acne-prone but also sensitive skin. Bio-effects: has a positive effect on acne-prone skin, reduces inflammation, regulates the functioning of the sebaceous glands and the amount of sebum secreted, the skin becomes tighter and enlarged pores become less visible.
4. Red photon light: Works at the level of the epidermis and dermis. Reduces and soothes all types of irritations. In addition, red light improves the elasticity and firmness of the skin. It delays the aging process and reduces wrinkles and discoloration. It is also great at reducing stretch marks, scars and cellulite.

5. **Green photon light:** Works at the epidermis level. Recommended for discoloration and sun spots. Penetrates into the basic layers of the epidermis, reduces expression lines, stimulates nutrients, and helps regenerate mature skin.

#### **Area of use:**

Derma Ice is equipped with 4 types of heads, each of which has a different purpose. Before using the massage heads, it is necessary to use a cosmetic that provides lubrication on the skin, e.g. cream, oil, serum or mask.

1. **Microcrystalline head:** can be used on any part of the face except the immediate eye area. It can also be used on areas of the body covered with hard, calloused skin, e.g. elbows or knees.
2. **Large round head:** can be used on any part of the face, except the eye area. Features: Powerful suction, suitable for larger surfaces (e.g. forehead). When using, move slowly over the skin and do not stop in one place. Swipe from the bottom to the top, moving slowly from the chin up.
3. **Small round head:** can be used on sensitive skin around the eyes, also on other parts of the face. Features: Weaker suction power. It can be used for areas where the skin is thin and delicate. Not suitable for strong massage.
4. **Oval head:** can be used on any part of the face and neck, except near the eyes. Features: Designed for lifting and firming, improves blood circulation in the skin and increases its elasticity.

#### **Service:**

1. **On/Off:** Press and hold the power on/off button for 3 seconds, the product is turned on, the display is on and the suction power displays 1 level. Then press the same button to activate subsequent suction levels.
2. **Power Adjustment:** Short press the mode switch button to adjust the suction power, each time you press it to increase one level and so on (1-2-3-4-5-0-1-2-3-4-5).
3. **HOT Mode:** Short press the mode switch button in the on state. The HOT mode light (sun icon) will turn on and switch to HOT mode automatically. Place the HOT surface on your face and move it over the entire face in circular motions to loosen the skin and open the pores.
4. **ICE Mode:** Press the mode switch button twice in the on state. The light (snowball icon) in ICE mode will turn on and switch to ICE mode automatically. Place the ICE surface on your face and move it over the entire area in circular motions to soothe the skin, close pores and regenerate it after microdermabrasion treatment.
5. **Photon light therapy:** Press three times will cut the mode switch button in the on state to switch to photon lights mode. After turning on this feature, move the device over your face in circular motions. From the inside out, from the bottom up, symmetrically left and right. Treat each area for approximately 30 seconds.
6. **Precautions:** Adjust the suction head to the 90-degree position and slowly move the suction head over your face. Don't hold in the same part for a long time. (Before removing the suction head from the skin, tilt it to the side to release the pressurized air. It can then be safely removed

without irritating the skin. Do not pull or move it directly to avoid pulling the skin and causing redness.)

**The device heats up after prolonged use of the ICE mode. Therefore, it is not recommended to use this mode for longer than 3-5 minutes. After use, turn off the product and wait 30 minutes before using it again.**

#### **Application:**

- discoloration
- scars
- stretch marks
- enlarged, very visible pores
- blackheads and acne (avoid inflammatory and purulent lesions)
- wrinkles
- skin devoid of healthy color with the first signs of photoaging

#### **How to use:**

##### **Microdermabrasion**

1. Cleanse and dry your face. This will relax the skin and open the pores. Then place the microcrystalline tip onto the large round cap of the device.
2. Turn on the device and place the tip on your face. Start with your chin and move the device from bottom to top. We recommend using blue photon light, which has antibacterial properties
3. Place the device on your nose and move from the inside to the outside, thoroughly cleaning each lobe of the nose.
4. Then place the device on your forehead and move it from the center to the outside in both directions.
5. In the last step, adjust the suction intensity and place the tip on your cheeks. Move the device from the inside to the outside. From the bottom up.
6. After completing the treatment, wash your face with a soothing liquid and use a moisturizer with a filter. We recommend using a moisturizing mask and using the ICE function together with photon light therapy.

##### **Lifting massage**

1. Use a cleansing liquid or peel-off mask and use the HOT mode for about 3-5 minutes to open the pores and relax the skin.
2. Place the appropriate massage tip on the device.
3. Turn on the photon light designed for the effects you want to achieve.
4. Point the tip at the skin at a 90-degree angle. You can use the level of suction that suits your skin's needs.
5. During the massage, move the device slowly and in the same direction each time - do not move the device in an oval shape. You can suck and release the skin where you need it. Do not use in areas of purulent acne or other skin inflammations.
6. Do not use the device on the same parts of the face for too long, as redness may appear on the skin. If this happens, please stop using for a while.
7. After the treatment, you can use the ICE mode to soothe the skin and close the pores.

**Tips:**

Do not wash the device with water. The Derma Ice device is not waterproof.

**Contraindications:**

- active purulent acne
- rosacea
- retinoid therapy
- acid treatments, lasers - at least 4 weeks interval before microdermabrasion
- psoriasis
- viral, bacterial and fungal skin diseases
- vascular skin
- tendency to scarring
- broken skin continuity (fresh wounds)
- skin cancer or suspicious moles in the treated area
- facial fillers - at least 4 weeks before microdermabrasion, otherwise the preparation may move.

**More information about the company and products [page 2 - f2]:**

- [www.beautifly.eu](http://www.beautifly.eu)
- [www.youtube.com/beautifly.eu](http://www.youtube.com/beautifly.eu)

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**[EN] Rules for handling electrical and electronic equipment:** Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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**[EN] Summary declaration of conformity:** The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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